

WHAT IS **PARKY**

Parky is a **Parkinson's Disease** remote patient monitor.



Parky received 510(k) clearance^{Rx} from the USFDA.





HOW IT **WORKS?**



HOW IT WORKS?



Clinician's decision making process is supported by objective data for disease management.



CLINICIAN

Reports are generated and shared with selected healthcare professional upon consent.

R_x

Patient downloads the app from iOS App Store for both iPhone and Apple Watch




Apple Watch monitors motor fluctuations on the wrist and transfers data to iPhone.



PARKY APP

Patients receive medication reminders both on iPhone and Apple Watch.


Weekly and monthly med compliance data is reported based on the reminder responses.



Stories enables Parky to update and share PD related content instantly.



Steps and total recording time is reported.



Patients can view their reports in real-time.



HOW TO **LEVERAGE PARKY?**



A COMPANION TOOL



Each Parky kit includes an iPhone, an Apple Watch and 10+ activation code cards to be distributed by the prescribing physician to the patients.





WHAT IS THE OUTPUT?

Parky generates patient symptom profiles, including the data types below:

- tremor
- dyskinesia
- steps
- recording time
- medication compliance

Reports are generated on a daily, weekly and monthly basis and automatically shared with the prescribing clinician.






Report Type : Daily

Patient: Jane Doe

20.06.2022



Legend

Tremor
Likelihood Percentage

Dyskinetic Symptom
Likelihood Percentage

Prescribed Medication Time

[Clinician Guide](#)

Daily Med Compliance

Medication Name	Compliance
Epamor 20mg	100%
Gepamor 20mg	100%

Med compliance is subject to patient's adherence.

Daily Total Recording Time

10 Hours

Recording time is subject to patient's hardware use.

Daily Total Steps

800 Steps

Steps count is subject to patient's hardware use.

Minute by Minute



FROM SNAPSHOT TO HIGH RESOLUTION

For research purposes, Parky can generate tailor-made reports based on collected raw data on tremor and dyskinesia episodes, mobility, and symptom correlation with medication compliance.



RESEARCH

Parky is build upon Movement Disorders API developed by Apple Inc.

SCIENCE TRANSLATIONAL MEDICINE | RESEARCH ARTICLE

BIOSENSORS

Smartwatch inertial sensors continuously monitor real-world motor fluctuations in Parkinson's disease

**Rob Powers¹, Maryam Etezadi-Amoli¹, Edith M. Arnold¹, Sara Kianian^{1,2},
Irida Mance¹, Maxsim Gibiansky¹, Dan Trietsch¹, Alexander Singh Alvarado¹,
James D. Kretlow¹, Todd M. Herrington^{3,4}, Salima Brillman⁵, Nengchun Huang⁶,
Peter T. Lin⁶, Hung A. Pham¹, Adeeti V. Ullal^{1*}**

Longitudinal, remote monitoring of motor symptoms in Parkinson's disease (PD) could enable more precise treatment decisions. We developed the Motor fluctuations Monitor for Parkinson's Disease (MM4PD), an ambulatory monitoring system that used smartwatch inertial sensors to continuously track fluctuations in resting tremor and dyskinesia. We designed and validated MM4PD in 343 participants with PD, including a longitudinal study of up to 6 months in a 225-subject cohort. MM4PD measurements correlated to clinical evaluations of tremor severity ($\rho = 0.80$) and mapped to expert ratings of dyskinesia presence ($P < 0.001$) during in-clinic tasks. MM4PD captured symptom changes in response to treatment that matched the clinician's expectations in 94% of evaluated subjects. In the remaining 6% of cases, symptom data from MM4PD identified opportunities to make improvements in pharmacologic strategy. These results demonstrate the promise of MM4PD as a tool to support patient-clinician communication, medication titration, and clinical trial design.

Apple conducted a clinical study for design and development of the algorithms for tremor and dyskinesia detection.

The study is known to be the largest in the field of PD wearables in terms of participant numbers and follow-up time.

List of clinical investigators:

- Salima Brillman, MD
*Parkinson's Disease and Movement Center of
Silicon Valley, Menlo Park, CA 94025, USA*
- Todd M. Herrington, MD
*Department of Neurology, Massachusetts
General Hospital, Boston, MA 02114, USA
Department of Neurology, Harvard Medical
School, Boston, MA 02115, USA*
- Nengchun Huang, MD
*Silicon Valley Parkinson's Center,
Los Gatos, CA 95032, USA*




BENEFITS

Less time
spent for disease
management
and decision
making



HCP


Better treatment titration



Objective Data




Increased medication adherence



PATIENT

Increased
quality of life



WE ARE h2o° therapeutics

h2o is a digital therapeutics developer with a focus on wearable technologies. By the use of AI, AR, and mobile technologies, we make real-time human data become a handy tool for disease management in certain therapeutic areas. Our core understanding of digital therapeutics lies in the power of AI, the robustness of their clinical health research, and high levels of user engagement.



FDA 510(K)
Cleared



Clinical Research
Phase II

covie

FDA BDD
Under Review



THANK YOU!

—

h2o° therapeutics

